

GE Healthcare

Care Plus Incubator Platform



For Fragile Lives

GE's Care Plus® incubator platform brings you a microenvironment that eases the newborn and helps you provide exquisite care. Each flexible, modular system is designed to make caregiving easier in a busy environment with varying procedures.

Comfortable and cozy

Active double-wall construction minimizes radiant heat loss and maintains a stable temperature – even with the door open for procedures.

Humidity systems provide variable humidity control to protect the newborn's skin.

Nurturing Care

Soft progressive alarms alert you without disturbing baby. Silenceable alarms beep softly as an initial warning, and change to a two-tone alarm after 30 seconds.

Accessibility and visibility.

Alarm readouts, control panels and reservoir are easily viewed. Caregiver and family have easy access to baby from all four sides through six portholes. Clear acrylic hood allows distortion-free viewing.



Building On Your Investment

Using a Dovetail Rail and shelf system as a framework, the Care Plus incubators can be configured to meet the demands of your care area. And are adaptable, allowing for upgrades for the future. A wide variety of options improve organization, save space and allow you to spend more time caring for newborn patients.

Standard Care Plus Configuration Includes

- Dovetail Rail system.
- Standard Humidifier provides adjustable humidity.
- Care Plus 3000: Three-position tilt for easy, precise adjustment to 0°, 3° and 6° at either end of the bed.
- Care Plus 4000: Continuous Tilt allows smooth, unlimited adjustment from 0° to 12°.

Care Plus Options:

- Servo Humidifier controls humidity from 40-85% at any incubator temperature.
- In-bed scale measures weight without moving the baby.
- Elevating base with foot-activated switch.
- Stationary base for mobility and storage at an economical price.
- In-frame drawer storage units
- Mounting accessories include a variety of posts and brackets for mounting ventilators, infusion pumps and IV lines.

GE has a passion for providing patients with the best in perinatal care. Our Care Plus incubators are designed with both the caregiver and the newborn in mind.



© 2006 General Electric Company – All rights reserved.

GE and GE Monogram are trademarks of General Electric Company.

Care Plus® is a trademark owned by Datex-Ohmeda, Inc. and are federally registered in the United States Patent Office.

GE Healthcare Finland Oy., doing business as GE Healthcare.

Healthcare Re-imagined

GE is dedicated to helping you transform healthcare delivery by driving critical breakthroughs in biology and technology. Our expertise in medical imaging and information technologies, medical diagnostics, patient monitoring systems, drug discovery, and biopharmaceutical manufacturing technologies is enabling healthcare professionals around the world discover new ways to predict, diagnose and treat disease earlier. We call this model of care „Early Health.“ The goal: to help clinicians detect disease earlier, access more information and intervene earlier with more targeted treatments, so they can help their patients live their lives to the fullest. Re-think, Re-discover, Re-invent, Re-imagine.

GE Healthcare
P.O. Box 900, FIN-00031 GE, Finland
Tel. +358 10 394 11 • Fax +358 9 146 3310

www.gehealthcare.com



GE imagination at work